

instavac II  
instavac II  
instavac II  
instavac II

Squire-Cogswell  
aeros®

# instavac® II



## Tabletop Aspirator

Ideal for:

Gastroenterology • O.R. backup • Hospital Code Carts • Surgicenters • O.R., E.R., and/or  
General Wall Vacuum Backup • Special Procedure Labs • Clinics • Physician's Offices

## instavac® II



**T**he Instavac II is the second generation of this cool running, table top aspirator. Its compact, lightweight design allows it to be mounted on crash carts or used for scoping or O.R. backup. Features include a fan-cooled rotary vane pump and rotary vacuum regulator. The Instavac II is widely used in gastroenterology, hospital code carts, surgicenters, O.R., E.R. and/or general wall vacuum backup, special procedure labs, clinics, and physician's offices.

### Features:

- ✓ Instant vacuum
- ✓ Fan cooled rotary carbon vane pump
- ✓ Vacuum & flow that surpass wall vacuum performance
- ✓ Safety overflow trap
- ✓ Perfect for G.I. labs, O.R. & doctors offices
- ✓ *3 year warranty*

### Available Accessories:

- ✓ Collection Canisters (Disposable & Glass)
- ✓ Aspirating Change Out Kits
- ✓ Filters & Tubing
- ✓ Additional Accessories

### Performance:

Vacuum Range:	0 to 550mm Hg
Free Air Flow:	38 LPM
Pump:	rotary carbon vane type

### Collection Device:

Canister:	disposable plastic with mechanical shutoff (adjustable bracket will accept any canister)
Capacity:	1200cc standard
Tubing:	8" with bacterial filter

### Controls:

Vacuum Regulator:	rotary type on panel
Vacuum Gauge:	calibrated in mm Hg

### Electrical Requirements:

AC:	120 volts, 60 Hz, 2A
	220 volts, 50 Hz, 2A

### Physical:

Weight:	15.5 lbs.
Dimensions:	8"H x 12"W x 9.5"D

