

NEW ... OptiVac® AC/DC PORTABLE SUCTION UNIT

Model G180



Tubing can be cut to desired length and attached to disposable hydrophobic/bacteria filter.

Model G180 with Disposable Hydrophobic/Bacteria Filter and optional Carrying Case

THE QUALITY OF GOMCO® WHERE YOU NEED IT MOST – The new G180 is the second generation OptiVac in the industry-leading Gomco line of portable suction equipment. Its compact design, ease of operation, rechargeable battery, AC power supply, and fully-enclosed unit are just a few of its advanced features.

FEATURES

- Run time operation on fully-charged battery for approximately 3 hours at maximum vacuum.
- 25-550 mm Hg vacuum range, coupled with greater than 30 lpm flow, provides powerful suction for faster aspiration.
- Fully chargeable in 8 hours.
- Will run on AC power even with low battery.
- Will run for approximately 1 hour after charged for just 2 hours.
- Long-life rechargeable battery eliminates the need for frequent replacement.
- Sealed water-resistant unit is easy to clean, and minimizes potential contamination or damage to the compressor and battery.
- Control panel at an angle for easy viewing.
- Quality-engineered thermoplastic provides impact resistance and durability.
- 3-year warranty on the unit, and 1 year on the battery.

MODEL	COLLECTION BOTTLE	MOTOR	COMPRESSOR	VACUUM RANGE	FLOW RATE	DIMENSIONS HxDxW	SHIPPING WEIGHT
G180	1500 ml Disposable Collection Canister	3.4 amps	Diaphragm	1-22 in. Hg 25-550 mm Hg	30 lpm minimum	9.4 x 7.5 x 16.8 in. 23.9 x 19.1 x 42.7 cm	11.4 lbs. 5 kg

ORDERING INFORMATION

G180	OptiVac Portable Suction Unit with Gauge and Regulator, ON & OFF Buttons, one plastic 1500 ml disposable collection canister with lid,	20-08-0003	Disposable Collection Canister, 1500 ml., stem inlet, 48/case
	one length of 15" clear flexible tubing (pump tube),	20-08-0016	Disposable Collection Canister, 1500 ml., stem inlet, 16/case
	one length of 72" clear flexible tubing (patient tube),	01-90-3928	Disposable Hydrophobic/Bacteria Filter, 1/4"x1/4" barb, 3/pkg.
	three hydrophobic/bacteria filters, AC power cord with hospital-grade plug	S610100	Suction Tubing Kit (one each 72" and 15" length, 1/4" ID, clear, PVC)
SU1060	Battery (12 volt DC, 5 amp)		
L200023	AC Power Cord		
L190-BAG	Carrying Case with shoulder strap		