

AARON 1250™-G OB/GYN TOTAL SYSTEM SOLUTION

We've combined our Aaron 1250 (see page 8) package with our superb smoke evacuation system (see page 12) and our latest mobile stand technology (see page 22) to produce a compact, mobile system for the office or hospital. As an added benefit, we included a loop electrode kit (5 of our most popular loops) and a smoke evacuator hose reducer fitting (see page 13) for attaching to a vaginal speculum to create the best GYN system on the market today.

STYLE

A1250G

STANDARD EQUIPMENT

- **A1250 electro-surgical generator**
 - + Three extra disposable electrodes
 - + One disposable pencil
 - + A1252C reusable grounding cord
 - + ESRE-1 five split grounding pads
- **SEGYN Smoke Shark system**
 - + SE01 Smoke Shark
 - + SF18 18-hour filter
 - + 786T 7/8" tube 6' long
 - + FSSE footswitch
 - + SERF reducer fitting
(attaches to vaginal speculum)
- **ESMS-C electro-surgical mobile stand**
 - + ESMS mobile stand, with instrument drawer
 - + ESMS-B bottom tray
- **ESLK electro-surgical loop kit**
 - + ES09 10 mm x 10 mm loop
 - + ES11 15 mm x 10 mm loop
 - + ES12 20 mm x 8 mm loop
 - + ES13 20 mm x 15 mm loop
 - + ES16 10 mm x 8 mm square

