


DISPOSABLE ELECTRODES FROM A NAME YOU TRUST!



Many Sizes Available

All Aaron® blade, ball, needle, loop, and square disposable electrodes use the standard 3/32" diameter stainless steel shafts. These electrodes, like all Aaron® electrodes, feature safety grip insulators, combining patient and user safety with easy insertion into and removal from the surgical pencil.

The loop and square electrodes feature tungsten wire for superior shape and wire integrity throughout the excision procedure, while minimizing thermal damage and yielding excellent histopathology specimens. Each molded loop and square electrode includes our patented tray to protect the loops from bending during shipping and handling.

Ball electrodes in 3mm and 5mm sizes are available for fulguration and desiccation, while needle electrodes are available for precise, pinpoint cutting and coagulation.

SuperCut™ Needles (ES 60-63) are made of Tungsten and keep their shape throughout the procedure.

The Aaron Arthroscopic electrodes are available in four tip styles and designed to fit all standard electro-surgical pencils.



LIST OF STYLES AND CONVERSION CHART

Loop Size Locator

	STRYKER/COLORADO	PREMIER	WALLACH	CONMED	VALLEY	ACMI/CI/RCON	UTAH	
SuperCur™ Tungsten Needles	ES60	3CM Super Fine Needle, 5/box	N101A					
	ES61	3CM Modified Super Fine Needle, 5/box	N103A		E1651			
	ES62	4.5CM Super Fine Needle, 5/box	N102A					
Needle Electrodes	ES63	4.5CM Modified Super Fine Needle, 5/box	N104A					
	ES02	Standard Needle, 25/box	COOPER	PREMIER	WALLACH	CONMED	VALLEY	
	ES03	Extended Needle, 25/box			909135	275-1502/9113	E1552	ACMI/CI/RCON
Dermal Loops	ES38	Modified Needle, 25/box				275-1504/9116	E1552-6	
	ES40	Extended Modified Needle, 25/box	N2512			138025		6162901
	ES22	1/4" Short Shaft Loop, 5/box						
	ES23	5/16" Short Shaft Loop, 5/box						
Dermal Tips	ES24	3/8" Short Shaft Loop, 5/box						
	ES25	1/2" Short Shaft Loop, 5/box						
	ES26	3/4" Short Shaft Loop, 5/box						
	A804	Sharp Dermal Tip, Non-Sterile, 100/box						
Blade Electrodes	A805	Sharp Dermal Tip, Sterile, 50/box						
	A806	Blunt Dermal Tip, Non-Sterile, 100/box						
	A807	Blunt Dermal Tip, Sterile, 50/box						
	ES01	Standard Blade, 50/box					275-1500/9111	E1551G
Ball Electrodes	ES04	Extended Blade, 25/box						E1551-6
	ES37	Modified Blade, 25/box						E1553
	ES39	Extended Modified Blade, 25/box						
	ES18	Anoled Blade, 5/box						
LLETZ Loop Electrodes	ES06	Extended 3mm Ball, 5/box	B0312			60-5181-102	E1563	6156901
	ES07	Extended 5mm Ball, 5/box	B0512			60-5181-103	E1564	6157901
	ES20	3/16" Ball, 25/box				275-1501/9112	E1550	
	ES21	5/32" Ball, 25/box						
Arthroscopic Electrodes	ES50	1/8" Ball, 25/box						
	ES08	5mm x 5mm Loop, 5/box						
	ES09	10mm x 10mm Loop, 5/box	R1010					6133901
	ES10	15mm x 8mm Loop, 5/box	R1507*					6134901
	ES11	15mm x 10mm Loop, 5/box						6135901
	ES12	20mm x 8mm Loop, 5/box	R2008			60-5180-102		6136901
	ES13	20mm x 15mm Loop, 5/box	R2015					6137901
	ES31	20mm x 20mm Loop, 5/box	R2020					6138901
	ES42	20mm x 12mm Loop, 5/box						
	ES43	15mm x 15mm Loop, 5/box						
	ES44	15mm x 5mm Loop, 5/box	R1515					
	ES45	13mm x 13mm Loop, 5/box	R1505					
	ES46	10mm x 5mm Loop, 5/box						
	ES47	25mm x 10mm Loop, 5/box						
	ES48	10mm x 3mm Loop, 5/box	R2510					
	ES49	20mm x 10mm Loop, 5/box						
ES51	13mm x 8mm Loop, 5/box	R2010			60-5180-104			
ES52	10mm x 8mm Loop, 5/box				60-5180-101*			
ES53	15mm x 12mm Loop, 5/box							
ES57	25mm x 8mm Loop, 5/box							
ES58	25mm x 22mm Loop, 5/box							
Square Electrodes	ES14	5mm x 5mm Square, 5/box						E1566
	ES15	10mm x 4mm Square, 5/box	S1004					6153901
	ES16	10mm x 8mm Square, 5/box						6154901
Arthroscopic Electrodes	ES41	5mm x 10mm Square, 5/box						6155901
	ES56	10mm x 10mm Square, 5/box						
	AR00	Electrode Hook 90, 5/box						
AR01	Electrode Hook 45, 5/box							
AR02	Electrode Acromioplasty, 5/box							
AR03	Electrode Meniscectomy, 5/box							

* shows 1mm difference in size

