

BARIATRIC SERIES

FEATURES

- Seat has additional depth and an extra center support for added strength
- Extra dimensions make it easier for patient and technician
- Adjustable height and depth flip arms
- Armrests lock in place when in use, then flips out of the way for easy exit
- Comfortable 2" padding throughout
- Easy-clean, antimicrobial, heavy weight, vinyl upholstery
- Neutral Grey powder-coated frame
- 2 knobs lock chrome armrest into position
- Adjustable foot levelers
- Chairs may be special ordered with straight armrests instead of flip arms (677 Option)
- Accommodates patients up to 700 lbs.



	Seat	Height
6000B-F	35" x 21"	20"
	Width	Depth
Overall	48"	32"
	Height Adjustment	
Arms	29"—34"	

Bariatric Blood Drawing Chair

- Accommodates larger patient or child with parent up to 700 lbs. combined weight
- Padded and upholstered seat, back and armrests
- 2 adjustable flip-arms

	Seat	Height
6000B-FD	35" x 21"	20"
	Width	Depth
Overall	58"	32"
	Height Adjustment	
Arms	29"—34"	

Bariatric Blood Drawing Chair with Drawer and Flip-Arm

- Accommodates larger patient or child with parent up to 700 lbs. combined weight
- Padded and upholstered seat, back and armrests
- 1 adjustable flip-arm
- Easy-clean, one-piece drawer

