

# MD SERIES

## FEATURES

- Contemporary all-welded steel frame
- Adjustable height and depth flip arm
- Armrest locks in place when in use, then flips out of the way for easy exit
- Comfortable 2" padding throughout
- Easy-clean, antimicrobial, heavy weight, vinyl upholstery
- White powder-coated frame
- 2 knobs lock chrome armrest into position
- Non-marking glides
- Adjustable arm width
- 300 lbs. capacity

	Seat	Height
<b>64950-F</b>	17" x 16"	20"
	Width	Depth
Overall	32"	26"
	Arm Width Adjustment	Arm Height Adjustment
Arm	19"—24"	28"—33"

### Blood Drawing Chair

- Padded and upholstered seat and armrests
- Adjustable arms
- White chair frame



Arm flips clear for easy seating and exit



Two knobs securely lock armrest into position for safety





	Seat	Height
<b>64929-F</b>	17" x 16"	20"
Overall	Width	Depth
	42"	26"
	Arm Width Adjustment	Arm Height Adjustment
Arm	19"—24"	28"—33"

**Blood Drawing Chair with Drawer**

- Padded and upholstered seat and armrests
- Adjustable arms
- White chair frame
- Easy-clean, one-piece plastic drawer



	Seat	Height
<b>64929-BF</b>	17" x 16"	20"
Overall	Width	Depth
	42"	26"
	Arm Width Adjustment	Arm Height Adjustment
Arm	19"—24"	28"—33"

**Blood Drawing Chair with Drawer**

- Padded and upholstered seat, back and armrests
- Adjustable arms
- White chair frame
- Easy-clean, one-piece plastic drawer



	Seat	Height
<b>64950-BF</b>	17" x 16"	20"
Overall	Width	Depth
	32"	26"
	Arm Width Adjustment	Arm Height Adjustment
Arm	19"—24"	28"—33"

**Blood Drawing Chair**

- Padded and upholstered seat, back and armrests
- Adjustable arms
- White chair frame