

**PUT YOUR CONFIDENCE IN...
AND BABIES ON...**

Detecto's Model 8450 Baby Scale



**Large, Easy-to-read, 1" LCD Display
Built-In Measuring Tape
Battery or AC Power Adapter
Memory and Lock-in Weight Features
Light Weight, One Piece Construction
Accurate Weights Displayed in lb/oz or kg**

DETECTO'S 8450 BABY SCALE...

Perfect for Hospitals, Physicians' Offices or, at Home

This extremely light weight scale can be used either in a permanent location or as a portable scale. The platform, base and padded carrying handles are incorporated into one complete unit making it easy to transport from place to place. Padded side rails and handle covers are slightly contoured ensuring baby will be cradled comfortably & securely while being weighed. Measure infant's length in inches or centimeters with built-in measuring tape.



One-piece,
lightweight
construction



Conveniently located,
easy-to-use display and keypad

The 8450 is extremely easy-to-use. Accurate weight is displayed in pounds/ounces or kilograms on the large, 1" LCD. Conversion is as simple as the touch of a button. Optional units, pounds only and ounces, are available via setup. The Memory feature allows operator to record weight. Lock-In weight display means baby's movements won't effect accurate weight readings.

Specifications	
Capacities:	20 kg x 0.01 kg/44 lb x 0.02 lb/704 oz x 0.5 oz
Weight Units:	lb/oz, kg, lb, oz
Power:	9V battery or 9V/100mA AC adapter (not included)
Display:	1" LCD display
Operation	
Temperature:	5° to 35°C/41° to 95°F
Keys:	On/Zero/Off, kg, lb, oz, lb/oz, hold
Measuring Tape:	Inch and centimeter graduations
Platform:	500 x 350 x 50 mm/20" x 14" x 2"
Overall	
Dimensions:	620 x 350 x 50 mm/24" x 14" x 2"
Net Weight:	3.5 kg
8450-CASE:	optional carrying case

We reserve the right to improve, enhance or modify features and specifications without prior notice.

