



Three Drawer Economy Treatment Cart #ETC-3



Specifications

PRODUCT REQUIREMENT

- Harloff's Economy Treatment Carts combine useful function with competitive price. These carts are ideal for budget sensitive end users, such as outpatient clinics, doctor offices, veterinary centers, dental offices and schools. Low height and small size make them easy to store under a counter or against a wall.

Construction

- Painted Steel Cabinet and Drawers

Standard Features

- Compact Size
- Keylocking system with two keys
- Roller Drawer Slides
- Removable Plastic Top
- Certified "Latex Free"
- Five Year Warranty
- 3-inch (7.6 cm) single-wheel stem caster set, full swivel – two with brake
- Hard-baked powder coat finish
- Cart requires caster installation after delivery

Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

Drawer configuration:

- 3 x 2.75-inch (7.0 cm)

Open Storage:

- 1 x 12 -inch (30.5 cm)

Assembled Cart Dimensions:

- 27.25" H x 23.38" W x 15.13" D (69.2 cm H x 59.4 cm W x 38.4 cm D)
- Cart Net Weight: 48 pounds (21.8 kg)
- Approx. Shipping Weight: 63 pounds (28.7 kg)
- **PLEASE NOTE THIS CART IS SHORTER THAN STANDARD COUNTER HEIGHTS**