



Shown with R130 standing sling.

Invacare® Reliant Stand-Up Lift

Model no. RPS350-1

Comfortable and secure, the Invacare Reliant Stand-Up lift is the answer to your patient transport needs. Reliant's ergonomic styling and easy maneuverability make it an ideal product for everyday patient handling. With its numerous safety features, the Reliant Stand-Up lift ensures caregiver security and resident peace of mind. Ideal for use with weight-bearing or fully dependent patients as well as those needing rehabilitation support, the Invacare Reliant Stand-Up lift will be a valuable addition to your safety-management system.

FEATURES

- Innovative adjustment features adapt to a wide variety of body sizes and shapes. That means comfort at any height, width and weight range up to 350 lb.
- Low friction casters significantly improve rollability and contributes to caregiver safety
- Multifunction slings allow for quick toileting, stand assist or full-support seated transfer
- New, adjustable pivoting leg support features a soft pad with lateral contour that secures legs in three positions for comfort at any height
- Wide, non-slip footplate provides stable base for transfer
- Under-bed clearance of 4.5" allows accessibility to most nursing home, hospital and home care bed styles
- Removable battery eliminates downtime
- Ergonomic features make it ideal for toilet transfers, as well as for dependent and rehabilitation purposes
- Features include fast, simple sling attachment that protects against accidental disengagement; wide push handle; good access for clothing removal; and comfortable slings with extended back support that prevents pressure in the axilla area
- Adjustable leg strap option secures leg against knee pad for optimum positioning on footplate
- Warranty includes three years on lift, one year on motor
- Manufactured in an ISO9001 facility



Lift Specifications


Dimensions	
Height	
Maximum at sling hookup	63.7"
Minimum at sling hookup	39.6"
Base height (clearance)	4.5"
Base width	
Open	36.8"
Closed	25.7"
Base length	35.5"
Casters	
Front	3"
Rear (locking)	5"
Weight capacity	350 lb.
Product weight	
Out of carton	115 lb.
Shipping	125 lb.

Sling Specifications


	R130	R134
Dimensions		
Width	44"	34"
Length	16"	42"
Product weight		
Out of carton	1.1 lb.	1.3 lb.
Shipping	1.6 lb.	1.8 lb.



Shown with R134 Transport Sling.



Model no. R134 Transport Sling



Model no. R130 Standing Sling

Complete instructions for lift usage and sling procedures are found in the Invacare Operating Manual. Only genuine Invacare slings and accessories should be used on Invacare lifts.



Electronics

- Removable battery box is conveniently located on mast; it features an emergency off button and recharge indicator*
- Quick-release battery eliminates downtime
- Anti-entrapment feature stops motion if boom meets resistance while lowering
- Heavy-duty hand pendant and cord
- Includes two batteries and a recharging station

*With audible recharge alarm.