

OEM & Special Order

Laboratory Incubators

- Variety of capacities available
- Simple keypad input
- Stainless steel interior measures 30x30x30cm
- Working temperature ambient +5° -70° C
- Adjustable shelves
- Double doors to allow sample viewing

Part#	Incubator Capacity
LW27I	27 Liters
LW50I	50 Liters
LW80I	80 Liters
LW150I	150 Liters
LW300I	300 Liters



Part# LW150I

Platelet Rotator

- Widely used in blood banks
- Smooth 360° rotation
- Platelet bags load easily in the stainless steel basket
- Built-in tachometer
- 4 capacities available

Part#	LWPR-2	LWPR-4	LWPR-6	LWPR-8
Number of Baskets	1	2	3	4
Capacity	12	24	36	48



Part# LWPR-6

Pipette Shaker

- Portable & compact to conserve space
- Capacity 6 or 8 pipettes
- 2500 rpm
- 15-minute timer
- Suction-cupped feet

Part#	Capacity
LW-PS6	6 pipette
LW-PS8	8 pipette



Part# LW-PS6